

PRESTAN PROFESSIONAL ADULT MEDIUM SKIN CPR-AED TRAINING MANIKIN (WITH CPR MONITOR) – USA

Demo:<https://www.youtube.com/watch?v=NwIED3DOFvk>

Brand : Prestan, USA



PRODUCT DESCRIPTION:

This manikin features a unique clamshell design that accommodates an easy-to-insert face shield lung bag, and Prestan's new CPR Rate Monitor that provides instant feedback to both students and instructors about the rate of chest compressions! Other features include a patented face/head tilt that simulates the way an actual victim's head would move, and helps teach students how to open the airway during CPR; full chest rise; and a clicker mechanism that helps students learn to use the right amount of force in order to compress the victim's chest to the correct depth. Portable and easy to set up at training sites, the manikin comes with 10 adult/child face shield lung bags.